



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT		Effective Date	04/12/2024
SITAGLIPTIN TABLET USP 50 MG			
Protocol No.	QCM/MVP/028	Page	1 of 14

PROTOCOL CUM REPORT FOR MLT (MICROBIAL LIMIT TEST) VALIDATION OF SITAGLIPTIN TABLET USP 50 MG



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT

Effective Date

04/12/2024

SITAGLIPTIN TABLET USP 50 MG

Protocol No.

QCM/MVP/028

Page

2 of 14

SR. NO.	INDEX	PAGE NO.
1.0	Protocol Approval	3
2.0	Objective	4
3.0	Purpose	4
4.0	Scope	4
5.0	Responsibility	4
6.0	Prerequisite	5
7.0	Procedure	5 to 8
8.0	Validation Requirement & Observation	8 to 12
9.0	Acceptance Criteria	12
10.0	Revalidation Criteria	13
11.0	Summary & conclusion	13
12.0	Term & Definition	14
13.0	Abbreviation	14-15
14.0	Reference	15
15.0	Annexures	15



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT		Effective Date	04/12/2024
SITAGLIPTIN TABLET USP 50 MG			
Protocol No.	QCM/MVP/028	Page	3 of 14

1.0 PROTOCOL APPROVAL:

Signing of this approval page of protocol indicates agreement with the qualification approach described in this document. If modification to the qualification approach becomes necessary, an addendum should be prepared and approved. The Protocol/Report cannot be used for execution unless approved by the following authorities.

PREPARED BY

NAME	DESIGNATION	SIGNATURE/DATE
Narain Mishra	Asst. manager	Rm 02/12/24

CHECKED BY

NAME	DESIGNATION	SIGNATURE/DATE
Sushmap. Nautiyal	Asst. Manager	Sushmap 02/12/24
Karmala Kumbh	Sr. Manager	Karmala 02/12/24

APPROVED BY

NAME	DESIGNATION	SIGNATURE/DATE
Karmala Kumbh	Dr. Manager	Karmala 03/12/24



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT

Effective Date

04/12/2024

SITAGLIPTIN TABLET USP 50 MG

Protocol No.

QCM/MVP/028

Page

4 of 14

2.0 OBJECTIVE

The objective of this study is to verify the method suitability of Microbial Enumeration Tests performed for determination of total aerobic microbial count (TAMC), total combined yeast and Mold count (TYMC), and free from specified microbial species in the pharmaceutical articles as per harmonized method.

3.0 PURPOSE

The purpose of this protocol is to give test method validation of microbial limit test method for **Sitagliptin tablet USP 50 mg**, acceptance criteria and summary.

4.0 SCOPE

This protocol is applicable for microbial enumeration tests and tests for specified microorganisms in **Sitagliptin tablet USP 50 mg** used in the manufacturing of Product at Next Wave (India), Paonta sahib.

5.0 RESPONSIBILITIES

The below mentioned departments are responsible for overall technical excellence, applicability in terms of Scheduling, execution, data compilation, completeness, accuracy report review and investigation (in case of failures) of simulation studies. The responsibilities of individual department of the validation group are defined as below

Department/ Function	Responsibility
Microbiology	Preparation of Protocol.
	Execution of validation activity as per the approved protocol.
	Documentation and completion of validation activity.
	Preparation of summary report.
	Review the activity and provide required technical support.
QA	To review the protocol and approved the protocol



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT		Effective Date	04/12/2024
SITAGLIPTIN TABLET USP 50 MG			
Protocol No.	QCM/MVP/028	Page	5 of 14

6.0 PREREQUISITE

- 6.1 Ensure the training is completed to all the executors involved in the study.
- 6.2 Ensure that the documents (COA, MSDS) are available (if any) for reference purpose.
- 6.3 Ensure that the Validation or Calibration status of the equipment's involved in the study.

7.0 PROCEDURE

7.1 Overview of the Procedure:

- 7.1.1 Important objectives of validation are determining.
 - 7.1.1.1 Whether the sample to be examined has any inherent antimicrobial properties.
 - 7.1.1.2 Whether the incubation and growth condition can recover microorganism that may be present to an acceptance level.
- 7.1.2 Neutralizing agents may be used to neutralize the activity of antimicrobial agents in products, see Table 1 for a list of potential neutralizing agents/methods. The appropriate neutralizing agent should be added preferably before sterilization of the media. Include a blank control with neutralizer and without product to demonstrate efficacy and absence of toxicity for microorganisms
- 7.1.4 Perform the method suitability before performing the Microbial enumeration tests and tests for specified microorganisms or simultaneously with the test for bioburden. The suitability of the test method shall be established preferably at the product development stage.
- 7.1.5 When method suitability test for Microbial enumeration tests and tests for specified microorganisms are performed simultaneously.
- 7.1.6 Summarize the method suitability data after completion of initial trials and conclusions shall be referred during routine testing and / or regulatory submission. Finding of trials/ runs and modifications done during study should be documented with appropriate rationale and conclusion.
- 7.1.7 Follow the test method given in below procedure.
- 7.1.8 Refer the Table 1 for selection of the technique and use of the neutralizers, if any, for recovery.



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT		Effective Date	04/12/2024
SITAGLIPTIN TABLET USP 50 MG			
Protocol No.	QCM/MVP/028	Page	6 of 14

Table 1

COMMON NEUTRALIZING AGENTS FOR INTERFERING SUBSTANCES

Interfering Substance	Potential Neutralizing Agents/ Method
Glutaraldehyde, mercurial	Sodium hydrogen sulfite (Sodium bisulfite)
Phenolics, alcohol, aldehydes, sorbate	Dilution
Aldehydes	Glycine
Quaternary ammonium compounds (QACs), Para hydroxybenzoates (parabens), bis to biguanides	Lecithin
QACs, iodine, parabens	Polysorbate
Mercurials	Thioglycolate
Mercurials, halogens, aldehydes	Thiosulfate
EDTA (edetate)	Mg or Ca ions

TYPICAL NEUTRALISING FLUID

Name of the ingredient of neutralizing	Concentration
Polysorbate 80	30 g
Lecithin (egg)	3 g
Histidine hydrochloride	1 g
Peptone (meat or casein)	1g
Sodium chloride	1g
Potassium dihydrogen phosphate	4.3 g
Disodium hydrogen phosphate dihydrate	7.2
Purified water /Water for injection	1000 ml

Sterilize by heating in an autoclave using validated cycle.

Note: If the solution has insufficient neutralizing capacity the concentration of polysorbate 80 /lecithin/any suitable ingredient may be altered (addition or deletion) based on product nature.

- The pre-suitability trials / runs should be elimination or neutralization of interfering factors there by establishing the required pH, diluent, quantity of neutralizer, number of rinsing / washing cycles to



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT		Effective Date	04/12/2024
SITAGLIPTIN TABLET USP 50 MG			
Protocol No.	QCM/MVP/028	Page	7 of 14

be employed for suitability study. In pre to suitability trails, neutralizer should be evaluated for its toxicity and efficacy.

7.2 Method suitability Study:

7.2.1 **Purpose of the study:** The method suitability study is performed to demonstrate that the Pharmacopeial article under test does not adversely affect the reliability of the test method with its inherent bacteriostatic and fungistatic properties.

7.3 Method suitability Approach

- 7.3.1 In this case pour plate method as applicable.
- 7.3.2 Select the appropriate Diluent and/ or Medium with or without neutralizer based on the product nature.
- 7.3.3 Neutralizing agents shall be used to neutralize the antimicrobial activity of the material being testing, wherever applicable.
- 7.3.4 These may be added to the chosen Diluent, Medium and/ or Rinsing fluid, preferably before sterilization.
- 7.3.5 General neutralizing agents for respective interfering agents (material under test) and the composition of typical neutralizing fluid, which can be potentially used for various materials, are provided in Table 1 for reference.
- 7.3.6 Based on the nature of the material, test conditions/ method and available knowledge, appropriate neutralizing agent/ fluid shall be selected and incorporated in to method suitability studies.
- 7.3.7 Perform all possible ways to get the satisfactory result and based on available data final conclusion shall be made.
- 7.3.8 For pour plate method, use culture suspension of equivalent concentration in such a way that the final solution shall have not more than 100 cfu/ml.

7.4 Materials / Equipment Requirement:

7.4.1 Autoclave, Laminar Air Flow, Incubators, Spatulas, Sterile Petri plates, Micropipettes and sterile tips, bottles/ Test tubes, Inoculation loop, Water bath, Balance, pH Meter, Vortex mixer.

7.5 Media and Reagents Requirements:

7.5.1 Soya bean casein digest agar, Soyabean casein digest medium, 0.1 % Normal Saline, Buffered Sodium chloride peptone solution ph 7.0 with 1% Polysorbate 80, Sabouraud dextrose agar, Cetrimide Agar, Mannitol Salt Agar.



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT		Effective Date	04/12/2024
SITAGLIPTIN TABLET USP 50 MG			
Protocol No.	QCM/MVP/028	Page	8 of 14

7.6 Preparation of Media:

- 7.6.1 Prepare microbiological media required for Microbial enumeration tests and tests for specified microorganisms as per current version of SOP No. QCM-030 "Procurement, receipt, storage, issuance, preparation, sterilization and usage of microbiology media".
- 7.6.2 Carry out the growth promotion test (GPT) for all the media used in the study. Pre to incubate the media tubes/plates as applicable. Alternatively, keep the negative and positive controls simultaneously along with the test sample wherever necessary.

7.7 Preparation of Test Strains:

- 7.7.1 Prepare the test strains as per the current version of SOP "Preparation of culture suspension and inoculum determination" or SOP No QCM-021 "Receipt, Storage, Handling and Pre-qualification of Microorganism.:
- 7.7.2 Use the test strains of not more than 5 passages removed from original master seed lot.
- 7.7.3 Prepare the following test strains for the method suitability study.

Name of the Strains	NCIM No.	Equivalent to ATCC No.
<i>Pseudomonas aeruginosa</i>	NCIM 2200	ATCC 9027
<i>Staphylococcus aureus</i>	NCIM 2079	ATCC 6538P
<i>Escherichia coli</i>	NCIM 2065	ATCC 8739
<i>Aspergillus brasiliensis</i>	NCIM 1196	ATCC16404
<i>Candida albicans</i>	NCIM 3471	ATCC 10231

8.0 VALIDATION REQUIREMENTS & OBSERVATIONS:

8.1 Pretreatment of Test Specimen:

- 8.1.1 Prepare the specimen to be tested, by treatment that is appropriate to its physical characteristics and does not alter the number and kind of microorganisms originally present, in order to obtain a solution or suspension of all or part of it in a form suitable for the test procedure to be carried out.
- 8.1.2 If necessary, adjust the pH of the solution to 6 to 8 by using sterile sodium hydroxide or hydrochloric acid.



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT

SITAGLIPTIN TABLET USP 50 MG

Effective Date

04/12/2024

Protocol No.

QCM/MVP/028

Page

9 of 14

8.1.4 In case of enumeration test performed on non to pharmacopeia material and/ or material potentially consist antimicrobial characteristics, the quantity of sample shall be less with similar or lower dilutions.

8.1.5 Various dilutions used for enumeration are described in below table 2.

Table 2

Dilution	Quantity of Test Specimen to be Transferred	Quantity of Diluent
1:10	10 gm or 10 ml of product	90ml of buffered sodium chloride peptone solution with or without Neutralizer
	1 gm or 1 ml of product	9ml of buffered sodium chloride peptone solution with or without Neutralizer
1:50	10 gm or 10 ml of product	490ml of buffered sodium chloride peptone solution with or without Neutralizer
	1 gm or 1 ml of product	49ml of buffered sodium chloride peptone solution with or without Neutralizer
1:100	10 gm or 10 ml of product	990ml of buffered sodium chloride peptone solution with or without Neutralizer
	1 gm or 1 ml of product	99ml of buffered sodium chloride peptone solution with or without Neutralizer

Note: In case the sample quantity tested is less than the compendia recommendation, then the results shall be represented accordingly as per the quantity tested.

8.2 Suitability of Counting Method (Pour Plate Method):

8.2.1 Use the following test strains for quantitative recovery. For recovery of each test strain, carry out the following controls. Dissolve 10 gm of Sitagliptin tablet USP 50 mg sample, unless otherwise specified, in sterile 90 ml buffered sodium chloride-peptone solution pH 7.0 with 1% Polysorbate 80.

8.2.2 **Product Control (PC):** Transfer 1 ml of pre-treated Sitagliptin tablet USP 50 mg test specimen dilution on to each of two sterile Petri dishes (used for TAMC & TYMC).

Add 20 to 25 ml of sterile Soyabean casein digest agar for bacteria and incubate at 30 to 35°C for 3 to 5 days and add 20 to 25ml of sterile Sabouroud dextrose agar for yeast and Mold plates and incubate at 20 to 25°C for 5 to 7 days.



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT		Effective Date	04/12/2024
SITAGLIPTIN TABLET USP 50 MG			
Protocol No.	QCM/MVP/028	Page	10 of 14

8.3 Product Positive Control (PPC):

8.3.1 Transfer 1ml of pre-treated Sitagliptin tablet USP 50 mg test specimen dilution on to each of Two sterile Petri dishes and add not more than 100 cfu of all challenging organisms as per table 3, Add 20 to 25 ml of sterile Soyabean casein digest agar for bacteria and incubate at 30 to 35°C for 3 to 5 days and add 20 to 25ml of sterile Sabouraud dextrose agar for yeast and Mold plates and incubate at 20 to 25°C for 5 to 7 days.

8.3.2 **Inoculum Control (IC):** Add the culture having NMT 100 CFU into the Petri plates in duplicate. Add 20 to 25 ml of sterile Soyabean casein digest agar for bacteria and incubate at 30 to 35°C for 5 days and add 20 to 25ml of sterile Sabouroud dextrose agar for yeast and Mold plates and incubate at 20 to 25°C for 5 to 7 days.

8.4 **Negative Control (NC):** Take 1 ml of sterile diluent on to sterile Petri plate Add 20 to 25 ml of sterile Soyabean casein digest agar for bacteria and incubate at 30 to 35°C for 5 days and add 20 to 25ml of sterile Sabouroud dextrose agar for yeast and Mold plates and incubate at 20 to 25°C for 5 to 7 days.

Table No. 3

8.5 Challenging Organism:

Name of the Strains	NCIM No.
<i>Pseudomonas aeruginosa</i>	NCIM 2200
<i>Staphylococcus aureus</i>	NCIM 2079
<i>Escherichia coli</i>	NCIM 2065
<i>Aspergillus brasiliensis</i>	NCIM 1196
<i>Candida albicans</i>	NCIM 3471

8.7 Interpretation of results:

8.7.1 Count the number of colonies obtained on each of the plates. Compare the counts of product positive control with the counts obtained in inoculum control.

8.7.2 Count the number of colonies and report the results as colony forming units (CFU) multiply with dilution factor per sample tested. In case of no growth, report the results as less than dilution factor or not Detected per sample tested.



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT		Effective Date	04/12/2024
SITAGLIPTIN TABLET USP 50 MG			
Protocol No.	QCM/MVP/028	Page	11 of 14

8.7.3 Record the validation observation in Annexure 1 and compare the recovery against the acceptance criteria.

8.7.4 If the recovery is not as per the acceptance criteria, modify the procedure and revalidate the modified procedure for recovery.

8.8 Suitability of Method for Specified Microorganisms:

Test for specified microorganisms shall be performed by direct inoculation method.

8.8.1 Test for *Escherichia coli*

8.8.1.1 **Product Control(PC):** Transfer 10 ml Pre-treated sample from point No.8.2.1 and add in 90ml of Soyabean casein digest medium (Represents 1g sample).

8.8.1.2 **Product Positive Control(PPC):** Transfer 10 ml Pre-treated sample from point No.8.2.1 and add in 90ml of Soyabean casein digest medium (Represents 1g sample). Then add *Escherichia coli* test strain containing not more than 100 CFU.

8.8.1.3 **Inoculum Control(IC):** Add *Escherichia coli* test strain containing not more than 100 CFU to 90ml of Soyabean casein digest medium.

8.8.1.4 **Negative Control(NC):** Transfer 10 mL of sterile chosen diluent to 90 mL of Soyabean casein digest medium as a Test Negative Control.

8.8.1.5 Incubate Product Positive Control and Positive Control containers at 30 to 35 ° C for 18 to 24 hours. Incubate Product Control and Negative controls containers at 30 to 35 ° C for 24 to 48 hours.

8.8.1.6 Shake the containers, transfer 1 mL of above culture medium to 5 different containers of 100ml of MacConkey broth and incubate at 42 to 44°C for 24 hours in case of product positive control and positive control (Inoculum control). In case of Test specimen and Negative control containers at 42 to 44°C for 48 hours.

8.8.1.7 After incubation, Subculture from each of the above containers on to each of the 5 MacConkey agar plates separately. Incubate Product Positive Control and inoculum control at 30 to 35°C for 72 hours. Incubate Test specimen and Negative Controls at 30 to 35°C for 72 hours. Upon incubation, if brick red colored colonies, surrounded by zone of precipitated bile are observed, this indicates possible presence of *Escherichia coli*.

If characteristic growth observed in the test specimen further identification tests to be performed. Record the method suitability observations in Annexure-1 and summarize.



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT		Effective Date	04/12/2024
SITAGLIPTIN TABLET USP 50 MG			
Protocol No.	QCM/MVP/028	Page	12 of 14

9.0 ACCEPTANCE CRITERIA:

- 9.1 The media used in the study shall comply with the growth promotion requirements.
- 9.1.1 All instruments / equipment's to be used in the study should be calibrated / validated for their defined set parameters.
- 9.1.2 Microbial Counts of each of the test strains obtained in Product Positive Control (in presence of the product) should be within factor of two from the counts obtained in Inoculum control.
- 9.1.3 There should be no growth of Microorganisms in Negative Controls.
- 9.1.4 The specified microorganisms shall be detected with the characteristic growth on selective media.
- 9.1.5 Test specimen results should comply with specification limits.

9.2 Recommendation in case of recovery failure:

- 9.2.1 This information shall serve to indicate that the sample is not likely to be contaminated with the given species of microorganisms.
- 9.2.2 In such events, the appropriate justification shall be included in the summary/ conclusions with sound scientific rationale, with reference to with available information from public domain.
- 9.2.3 However, monitoring should be continued in order to establish the spectrum of inhibition and bactericidal / fungicidal activity of the article at a dilution where nearest recovery is obtained.
- 9.2.4 Failure of the organism to grow in the relevant medium invalidates that portion of the examination and necessitates a modification of the procedure by.
- 9.2.5 An increase in the volume of diluents or culture medium.
- 9.2.6 Incorporate of a specific or general neutralizing agent into the diluent.
- 9.2.7 Membrane filtration; or a combination of the above measures.
- 9.2.8 Then perform the test with the highest dilution factor compatible with microbial growth and the specific acceptance criteria.

10.0 REVALIDATION CRITERIA:

- 10.1 Revalidation shall be performed
- 10.1.1 When there is a change in the experimental conditions.

11.0 SUMMARY & CONCLUSION

- 11.1 Results shall be documented in the test data sheets provided as Annexure-1 to the protocol.
- 11.2 Based on the observations recorded, evaluation of the results shall be carried out.



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT

SITAGLIPTIN TABLET USP 50 MG

Effective Date

04/12/2024

Protocol No.

QCM/MVP/028

Page

13 of 14

11.3 On the basis of the results and evaluation, Conclusion shall be updated in record of observation. The summary report shall signify a conclusion to confirm the successful qualification of the procedure (or) revalidation (in case of failure) is required. Any modification in the method required before initiating revalidation shall be included in the report. Method adapted for sample preparation to overcome any interference shall be indicated in summary so that the same can be followed in routine analysis. After completion of validation, it shall be concluded that the average number of CFU recovered from the challenged test preparation must not differ by a factor greater than 2 from the calculated value for standardized inoculum.

12.0 TERMS AND DEFINITIONS:

Term	Definition
Aerobic organisms	Microorganisms which utilize oxygen as the final electron acceptor during metabolism; microorganism which will grow only in the presence of oxygen.
Bioburden	Total number of viable microorganisms on or in a health care product prior to sterilization.
Colony forming unit	Visible outcome of growth of microorganisms arising from a single or multiple cell.
Environmental flora	Microorganisms associated with a processing environment.
Growth promotion test	Test performed to demonstrate that media will support microbial growth.
Positive sample	Test sample which exhibits detectable microbial growth after incubation

13.0 ABBREVIATIONS:

Abbreviation	Description
SCDM	Soya bean casein digest medium
QA	Quality Assurance
SOP	Standard Operating Procedure
QA	Quality Assurance
HVAC	Heating Ventilation and Air Conditioning System
GPT	Growth Promotion Test
Ltr.	Liters
±	Plus, or minus
QC	Quality Control



NEXT WAVE (INDIA)

RAMPUR GHAT ROAD, PAONTA SAHIB, DISTT. SIRMOUR (H.P)

PROTOCOL CUM REPORT FOR MLT VALIDATION

NAME OF PRODUCT

SITAGLIPTIN TABLET USP 50 MG

Effective Date

04/12/2024

Protocol No.

QCM/MVP/028

Page

14 of 14

Abbreviation	Description
°C	Degree Centigrade
GTP	General Test Procedure
SS	Stainless Steel
mL	Milliliter
mm	Millimeter
% W/V	Percent Weight by volume
Gm	Gram
CFU	Colony Forming Unit
TYMC	Total Yeast and Mold Count
TAMC	Total Aerobic Microbial Count
USP	United States of Pharmacopeia
EP	European Pharmacopoeia
MET	Microbial Enumeration Test
ROB	Record of Observation

14.0 REFERENCES

- 14.1 EP to European Pharmacopoeia –General Chapter <2.6.1.2>
- 14.2 USP to United States of Pharmacopoeia to General Chapter <61, 62>, <1225> & <1227>
- 14.3 USP to Microbial Examination of non to sterile products: Acceptance criteria for pharmaceutical preparations and substances for pharmaceutical use <1111>

15.0 ANNEXURES

S. NO.	NAME OF ANNEXURE	ANNEXURE NO.
1	Report for MLT (microbial limit test) validation of Sitagliptin tablet USP 50 mg.	Annexure-1

END OF THE DOCUMENT